



ATMOS MedDoc 2.6

The new Standard in Video recording and Documentation



▶ Standard AVI uncompressed High quality videos and pictures.



▶ Easy management of patients incl. findings and diagnoses in unlimited database



▶ Video Editing capabilities.



▶ Easy to use export options to Word documents, USB drives, CD's or Email.



▶ Works on every standard Laptop/PC with Firewire. No special equipment needed.

Visualization and Documentation of Images, Diagnoses and Findings in The ENT Practice

ATMOS MedDoc 2.6 can be used with any voice, speech and swallowing evaluation system, whether it is mobile or stationary. MedDoc can simplify your practice by keeping your work organized and saving you time. This software provides perfect image quality by recording uncompressed videos that are easily created, marked, edited, and exported. Additionally, **MedDoc 2.6** provides a library of common ENT diagnoses` and findings that can be easily administered to a patient database. MedDoc software integrates easily with existing components, enabling you to quickly begin recording and exporting pictures and videos with MedDoc's advanced new features.



Manage your patients with the **MedDoc Database**


Time Saving

Most documentation systems have poor image quality and produce piles of unsorted print outs, video tapes, or CDs. With **ATMOS MedDoc 2.6**, you can create videos and pictures quickly and easily that are automatically archived to the system. With the ability to delete, edit, or add a video, you can easily find or create the documents you are looking for which eliminates clutter, keeps your work organized, and saves time.



Record a video and create a report in less than 5 minutes.





Ask for a free 30 day trial version of MedDoc 2.6

User Friendly

ATMOS MedDoc is very easy to follow and operate. It provides the option of recording by using either a foot pedal, mouse, or keypad which is similar to any video recorder. Videos are stored automatically and can be easily accessed using the patient database which can sort by either name or ID-number. You can also choose from a library of findings and diagnoses. **MedDoc 2.6** comes with a print out user manual containing step by step instructions with screenshots in addition to a video user manual which provides click by click directions.

Digital Archiving

The digital archiving system enables you to save what is important, and throw the rest away. If you are not satisfied with a recording, you can quickly record a new one and delete the first. **MedDoc 2.6** also has the feature of editing the video so you can cut out only the important part and discard the remainder. This allows you to save storage space and to avoid extra clutter of unnecessary information so you can find your edited pictures and videos easily.

Easy Exporting

With **Med Doc 2.6**, it is easy to export photos, films, and diagnoses/ findings to standard data formats, to internal or external storage. You can easily save a picture file, video file, or word document onto your computer, and then email it. You also have the option of burning a video to a DC or DVD, or save any type of file on a USB drive or external hard drive. This is ideal for exchange with patients, colleagues, or for giving presentations and lectures. Because **MedDoc 2.6** uses only uncompressed standard AVI Video, you can play the videos on any computer without additional software.

| Technical Specifications | |
|--|---|
| Minimum System Requirements | <ul style="list-style-type: none"> ▶ Intel Pentium 4 with 3 GHz or Centrino 1.7 GHz ▶ 512 MB RAM (1 GB or greater recommended) ▶ IEEE 1394 port (Firewire) ▶ 160 GB Hard Drive* ▶ CD/DVD burner ▶ Display 1042 x 768 (32 bit) ▶ Microsoft Windows® XP SP2 or Windows® Vista ▶ Direct X 9.0b |
| Disc Space for a 1 minute Video sequence | circa 120 MB |
| Videoformat | Uncompressed AVI |

* 750 GB or greater external hard drive for FEES procedures recommended

ATMOS, Inc.
 3717 Huckleberry Road
 Allentown, PA 18104
 Toll free: 1-866-862-8667
 Fax: (610) 351-6827
 sales@atmosmed.com
 www.atmosmed.us

Subject to technical alterations
 MEDDOC.V1.R0.081908.US